

Anti-PAI-1 antibody [MA-33H1F7] (STJ16100781)

STJ16100781

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-PAI-1 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Mouse
Reactivity	Human, Mouse, Rabbit, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MA-33H1F7
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	5054
Gene Symbol	SERPINE1
Uniprot ID	PAI1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	