

Anti-Galectin-3 antibody [B2C10] (STJ16100782)

STJ16100782

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Galectin-3 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, FC, ELISA, IHC-P, WB
Host/Source	Mouse
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B2C10
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50 IHC 1:50 ELISA 1:50 WB 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3958
Gene Symbol	LGALS3
Uniprot ID	LEG3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	