

Anti-SR-BI antibody [3D12] (STJ16100785)

STJ16100785

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-SR-BI is suitable for use in Immunohistochemistry and Flow Cytometry research applications. |
| Applications | IHC-F, FC |
| Host/Source | Mouse |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 3D12 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.1% BSA and 0.02% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|-------------------------|-------------|
| Gene ID | 949 |
| Gene Symbol | SCARB1 |
| Uniprot ID | SCRB1_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |