

Anti-Lactoferricin Bovine antibody [5F12.1.2] (STJ16100786)

STJ16100786

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lactoferricin Bovine is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Mouse
Reactivity	Bovine

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5F12.1.2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	