

Anti-HSV1 antibody [T96] (STJ16100788)

STJ16100788

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Herpes simplex virus type 1 is suitable for use in Western Blot, Immunocytochemistry and ELISA research applications.
Applications	WB, ICC, ELISA
Host/Source	Mouse
Reactivity	HSV-1

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	T96
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence