

Anti-Lactoferricin antibody [a-bC-lobe] (STJ16100790)

STJ16100790

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Bovine Lactoferricin is suitable for use in ELISA, Immunohistochemistry and Western Blot research applications.
Applications	ELISA, IHC-P, WB
Host/Source	Mouse
Reactivity	Bovine, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	a-bC-lobe
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence