

Anti-CNF1 antibody [NG8] (STJ16100792)

STJ16100792

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CNF1 is suitable for use in ELISA and Western Blot research applications. |
| Applications | ELISA, WB |
| Host/Source | Mouse |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | NG8 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.1% BSA and 0.02% sodium azide |
| Isotype | IgG2a |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| |
|---------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen Sequence |