

Anti-Cellular Fibronectin antibody [DH1] (STJ16100796)

STJ16100796

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cellular Fibronectin is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-F, WB, ELISA
Host/Source	Mouse
Reactivity	Human, Chicken, Rat, Rabbit, Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DH1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence