

Anti-Tenascin-C antibody [EB2] (STJ16100797)

STJ16100797

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tenascin-C is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-F, WB, ELISA
Host/Source	Mouse
Reactivity	Human, Rabbit, Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EB2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3371
Gene Symbol	TNC
Uniprot ID	TENA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	