

Anti-Tenascin-C antibody [DB7] (STJ16100798)

STJ16100798

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tenascin-C is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P, IHC-F, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DB7
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3371
Gene Symbol	TNC
Uniprot ID	TENA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	