

Anti-Plasma Fibronectin antibody [BF12] (STJ16100801)

STJ16100801

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Plasma Fibronectin is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BF12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	