

Anti-Vitronectin antibody [BE10] (STJ16100802)

STJ16100802

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Vitronectin is suitable for use in Immunohistochemistry and Western Blot research applications. |
| Applications | IHC-F, WB |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|-------------------------------------|
| Clonality | Monoclonal |
| Clone ID | BE10 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 1% BSA & 0.1% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 7448 |
| Gene Symbol | VTN |
| Uniprot ID | VTNC_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |