

Anti-Beta1-Integrin antibody [DF5] (STJ16100803)

STJ16100803

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Beta1-Integrin is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-P, IHC-F, WB, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DF5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P IHC-F WB (1:1000) ELISA
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	