

Anti-Beta1-Integrin antibody [DF7] (STJ16100804)

STJ16100804

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Beta1-Integrin is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-P, IHC-F, WB, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DF7
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P IHC-F WB (1:1000) ELISA
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence