

Anti-Epidermal Growth Factor antibody [F5] (STJ16100806)

STJ16100806

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Epidermal Growth Factor is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-P, IHC-F, WB, ELISA
Host/Source	Mouse
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	F5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P (1:20) IHC-F (1:20) WB ELISA
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	