

Anti-Intra-acrosomal protein antibody [Hs-8] (STJ16100809)

STJ16100809

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Intra-acrosomal protein is suitable for use in Western Blot, Immunocytochemistry and Flow Cytometry research applications.
Applications	WB, ICC, FC
Host/Source	Mouse
Reactivity	Pig, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Hs-8
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence