

Anti-Intra-acrosomal protein antibody [Hs-14] (STJ16100810)

STJ16100810

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Intra-acrosomal protein is suitable for use in Flow Cytometry, Western Blot and Immunocytochemistry research applications.
Applications	FC, WB, ICC
Host/Source	Mouse
Reactivity	Mouse, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Hs-14
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence