

Anti-Alpha-Fodrin antibody [AA6] (STJ16100815)

STJ16100815

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alpha-Fodrin is suitable for use in Immunohistochemistry, Western Blot and ELISA research applications.
Applications	IHC-F, WB, ELISA
Host/Source	Mouse
Reactivity	Human, Rat, Chicken, Rabbit, Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AA6
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F WB (1:1000) ELISA
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence