

Anti-Cytokeratin 8 antibody [2A4] (STJ16100817)

STJ16100817

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cytokeratin 8 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2A4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F WB (1:1000)
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3856
Gene Symbol	KRT8
Uniprot ID	K2C8_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	