

## Anti-Cytokeratin 7 antibody [4F5] (STJ16100818)

STJ16100818

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Cytokeratin 7 is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC-F, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	4F5
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-F WB (1:1000)
<b>Formulation</b>	PBS with 1% BSA & 0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3855</a>
<b>Gene Symbol</b>	<a href="#">KRT7</a>
<b>Uniprot ID</b>	<a href="#">K2C7_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	