

Anti-Vimentin antibody [65E] (STJ16100819)

STJ16100819

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Vimentin is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	65E
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7431
Gene Symbol	VIM
Uniprot ID	VIME_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	