

Anti-CD9 antibody [MM2/57] (STJ16100822)

STJ16100822

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD9 is suitable for use in Western Blot, Flow Cytometry and Immunohistochemistry research applications.
Applications	WB, FC, IHC-F
Host/Source	Mouse
Reactivity	Human, Llama, Mink, Mustelid, Ferret, Dog, Horse, Bovine, Cat, Rhesus Monkey, Pig, Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MM2/57
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	928
Gene Symbol	CD9
Uniprot ID	CD9_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	