

Anti-CD19 antibody [LE-CD19] (STJ16100823)

STJ16100823

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD19 is suitable for use in Western Blot, Immunohistochemistry, Immunoprecipitation, Flow Cytometry and ELISA research applications.
Applications	WB, IHC-P, IP, FC, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LE-CD19
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	930
Gene Symbol	CD19
Uniprot ID	CD19_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	