

Anti-Alkaline Phosphatase antibody [V17.1] (STJ16100829)

STJ16100829

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alkaline Phosphatase is suitable for use in ELISA and Immunohistochemistry research applications.
Applications	ELISA, IHC-F
Host/Source	Mouse
Reactivity	Human, Cow

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	V17.1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

