

Anti-GnRH Receptor antibody [F1G4] (STJ16100836)

STJ16100836

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-GnRH Receptor is suitable for use in ELISA, Immunohistochemistry and Western Blot research applications.
Applications	ELISA, IHC-F, IHC-P, WB
Host/Source	Mouse
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	F1G4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2798
Gene Symbol	GNRHR
Uniprot ID	GNRHR_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

