

Anti-SNAP-25 26-27kDa antibody [SP12] (STJ16100840)

STJ16100840

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SNAP-25 26-27kDa is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB, IHC-P, ELISA
Host/Source	Mouse
Reactivity	Human, Rat, Pig, Hamster, Gerbil

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SP12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence