

## Anti-CD55 antibody [67] (STJ16100845)

STJ16100845

### GENERAL INFORMATION


|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Mouse monoclonal antibody anti-CD55 is suitable for use in Western Blot, Flow Cytometry and Immunohistochemistry research applications. |
| <b>Applications</b>      | WB, FC, IHC-F   |
| <b>Host/Source</b>       | Mouse   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                            |                              |
|----------------------------|------------------------------|
| <b>Clonality</b>           | Monoclonal                   |
| <b>Clone ID</b>            | 67                           |
| <b>Concentration</b>       |                              |
| <b>Conjugation</b>         | Unconjugated                 |
| <b>Purification</b>        |                              |
| <b>Dilution Range</b>      |                              |
| <b>Formulation</b>         | PBS with 0.09% sodium azide  |
| <b>Isotype</b>             | IgG1                         |
| <b>Storage Instruction</b> | Store at 2-8°C upon receipt. |

### TARGET INFORMATION

|                           |           |
|---------------------------|-----------|
| <b>Gene ID</b>            | 1604      |
| <b>Gene Symbol</b>        | CD55      |
| <b>Uniprot ID</b>         | DAF_HUMAN |
| <b>Immunogen</b>          |           |
| <b>Immunogen Region</b>   |           |
| <b>Specificity</b>        |           |
| <b>Immunogen Sequence</b> |           |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081