

Anti-Mcm2 antibody [CRCT2.1] (STJ16100851)

STJ16100851

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Mcm2 is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB, IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CRCT2.1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4171
Gene Symbol	MCM2
Uniprot ID	MCM2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	