

Anti-Mcm5 antibody [CRCT5.1] (STJ16100852)

STJ16100852

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Mcm5 is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB, IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CRCT5.1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4174
Gene Symbol	MCM5
Uniprot ID	MCM5_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	