

## Anti-Actin Alpha-muscle antibody [HHF35] (STJ16100854)

STJ16100854

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Actin Alpha-muscle is suitable for use in EM, ELISA, Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	EM, ELISA, ICC, IHC-F, IHC-P, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Chicken, Monkey, Rabbit, Rat, Swine, Zebrafish

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	HHF35
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen Sequence</b>



---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081