

Anti-TGF Beta antibody [TB21] (STJ16100866)

STJ16100866

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TGF Beta is suitable for use in Western Blot, Immunohistochemistry, ELISA, Flow Cytometry and Immunohistochemistry research applications.
Applications	WB, IHC-P, ELISA, FC, IHC-F
Host/Source	Mouse
Reactivity	Human, Mink, Sheep, Rabbit, Rat, Mustelid, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TB21
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence