

Anti-SKP2 antibody [SKP2-8D9] (STJ16100869)

STJ16100869

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SKP2 is suitable for use in Western Blot, Immunohistochemistry, Immunocytochemistry, Immunofluorescence, ELISA and Immunoprecipitation research applications.
Applications	WB, IHC, ICC, IF, ELISA, IP
Host/Source	Mouse
Reactivity	Human, Primates

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SKP2-8D9
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB IHC ICC IF (1:250) ELISA IP
Formulation	PBS with 0.1% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6502
Gene Symbol	SKP2
Uniprot ID	SKP2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	