

Anti-Porcine Acrosin antibody [ACR.2] (STJ16100885)

STJ16100885

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Porcine Acrosin is suitable for use in Flow Cytometry, Immunocytochemistry and Western Blot research applications.
Applications	FC, ICC, WB
Host/Source	Mouse
Reactivity	Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ACR.2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

