

Anti-Neural Cell Adhesion Molecule antibody [RNL-1] (STJ16100890)

STJ16100890

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Neural Cell Adhesion Molecule is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RNL-1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F (1:10-1:20) WB
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4684
Gene Symbol	NCAM1
Uniprot ID	NCAM1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	