

Anti-LRP/MVP antibody [LMR5] (STJ16100893)

STJ16100893

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-LRP/MVP is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F, IHC-P, FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LMR5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	supernatant with 1% BSA and 0.1% sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	9961
Gene Symbol	MVP
Uniprot ID	MVP_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	