

Anti-Nitrotyrosine antibody [HM11] (STJ16100894)

STJ16100894

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Nitrotyrosine is suitable for use in Immunohistochemistry, ELISA, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, ELISA, IHC-P, WB
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HM11
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence