

Anti-MRP3 antibody [M3II-9] (STJ16100896)

STJ16100896

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MRP3 is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	ICC, IHC-F, WB, FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M3II-9
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	8714
Gene Symbol	ABCC3
Uniprot ID	MRP3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	