

Anti-MRP2 antibody [M2II-12] (STJ16100898)

STJ16100898

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MRP2 is suitable for use in Immunocytochemistry, Immunohistochemistry, Flow Cytometry and Western Blot research applications.
Applications	ICC, IHC-F, FC, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M2II-12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ICC (1:20-1:50) IHC-F (1:20) FC WB (1:50)
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	1244
Gene Symbol	ABCC2
Uniprot ID	MRP2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	