

Anti-MVP/LRP antibody [1027] (STJ16100911)

STJ16100911

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MVP/LRP is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-F, IF, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1027
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgM
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	9961
Gene Symbol	MVP
Uniprot ID	MVP_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

