

Anti-Breast Cancer Resistance Protein antibody [BXP-53] (STJ16100917)

STJ16100917

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Breast Cancer Resistance Protein is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	ICC, IHC-F, WB, FC
Host/Source	Rat
Reactivity	Mouse, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BXP-53
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ICC 1:20-1:50 IHC 1:20 WB 1:20-1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	9429
Gene Symbol	ABCG2
Uniprot ID	ABCG2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	