

Anti-CD13 antibody [WM15] (STJ16100918)

STJ16100918

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD13 is suitable for use in Flow Cytometry, Immunoprecipitation and Immunohistochemistry research applications.
Applications	FC, IP, IHC-F
Host/Source	Mouse
Reactivity	Human, Primates

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	WM15
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	290
Gene Symbol	ANPEP
Uniprot ID	AMPN_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	