

Anti-Endothelial Cell Protein Receptor antibody [LMR-42] (STJ16100919)

STJ16100919

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Endothelial Cell Protein Receptor is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F, FC
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LMR-42
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F (1:20) FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence