

Anti-MRP2 antibody [M2III-5] (STJ16100920)

STJ16100920

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MRP2 is suitable for use in Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
Applications	ICC, IHC-F, IHC-P, WB
Host/Source	Mouse
Reactivity	Rat, Mouse, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M2III-5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ICC IHC-F IHC-P (1:20) WB
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	1244
Gene Symbol	ABCC2
Uniprot ID	MRP2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	