

Anti-CD120b antibody [80M2] (STJ16100924)

STJ16100924

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD120b is suitable for use in Flow Cytometry, ELISA, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, ELISA, IF, IP
Host/Source	Mouse
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	80M2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	