

Anti-TNF-Alpha antibody [4H31] (STJ16100927)

STJ16100927

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TNF-Alpha is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F, FC, ELISA, IP, WB
Host/Source	Mouse
Reactivity	Human, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4H31
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7124
Gene Symbol	TNF
Uniprot ID	TNFA_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	