

Anti-CD120a antibody [MR1-2] (STJ16100928)

STJ16100928

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD120a is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MR1-2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	