

Anti-TNF-RI antibody [HM104] (STJ16100930)

STJ16100930

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-TNF-RI is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA and Immunoprecipitation research applications.
Applications	IHC-F, FC, ELISA, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HM104
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7132
Gene Symbol	TNFRSF1A
Uniprot ID	TNR1A_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	