

Anti-TNF-RII antibody [HM102] (STJ16100931)

STJ16100931

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-TNF-RII is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F, FC, ELISA, IP, WB
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HM102
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7133
Gene Symbol	TNFRSF1B
Uniprot ID	TNR1B_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	