

Anti-TNF-Alpha antibody [T1] (STJ16100934)

STJ16100934

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TNF-Alpha is suitable for use in Flow Cytometry, ELISA and Western Blot research applications.
Applications	FC, ELISA, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	T1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7124
Gene Symbol	TNF
Uniprot ID	TNFA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	