

Anti-MRP4 antibody [M4I-80] (STJ16100936)

STJ16100936

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP4 is suitable for use in Immunocytochemistry, Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	ICC, IHC-F, WB, FC
Host/Source	Rat
Reactivity	Human, Mouse, Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M4I-80
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ICC IHC-F (1:20) WB (1:50) FC
Formulation	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	10257
Gene Symbol	ABCC4
Uniprot ID	MRP4_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	